

## Patient protection eyewear "ALLROUND", color silver

### Product information:

Item no.: 600-ALLROUND-20  
Product weight: 0 g



### Certified protection specifications for EN 207

#### Wavelengths

- 180-315
- >315-1400
- >1400-11000

#### Protection levels

D LB10 + IR LB5 (OD10+)  
D LB7 + I LB9 + R LB8 + M LB10 (OD10+)  
DI LB5 (OD5+)  
PF S CE

## Properties:

Patient goggles for almost any laser/wavelengths, made of stainless steel with opaque, non-reflective surface material for maximum safety. Very high wearing comfort due to variably adjustable headband.

The additional rim rubbers ensure a firm and comfortable hold.

- Application: Patient goggles (eyeball goggles) for almost any laser/wavelengths
- Material: stainless steel, non-reflecting surface (opaque)
- Color: silver
- with turquoise silicon rim for improved wearing comfort, plus 1 spare pair enclosed

Included in delivery: ALLROUND-20 patient goggles, eye rings (4x), length-adjustable eyewear cord in black, hard plastic case, instructions for use.