

Eyecaps metal incl. suction cup

Product information:

Item no.: 600-AUGK-M-XS

Product weight: 0 g

Application for all laser wavelengths Material:
stainless steel, opaque Dimensions: approx. 22.0 x
19.5 mm A pair with a rubber sucker Eye caps
autoclavable



Properties:

Our **eye caps out of metal** convince with a high wearing comfort and are available in different sizes.

Application for all laser wavelengths

Material: stainless steel, opaque

Dimensions: approx. 22.0 x 19.5 mm

A pair with a rubber sucker

Eye caps autoclavable

EYE CAPS ARE NOT DISCOUNTABLE!