



IPL protection eyewear ONTOR F-4

Product information:

Item no.:	600-K234-ONTO-54
Filter color:	grey
Frame color:	black/white
Filter material:	polycarbonate, scratch resistant
Product weight:	44 g



Alle Angaben ohne Gewähr.
All information is provided without guarantee.

IPL safety eyewear ONTOR F4 (suitable also for spectacle wearers) offer safety for operator and patient during IPL treatments.

Properties:

Our **IPL safety goggles ONTOR** are lightweight, versatile **safety goggles** made of a plastic frame, which is perfectly suitable as over-glasses for spectacle wearers. The frame offers a very large field of vision and wide side coverage, with excellent wearing comfort due to flexible temples adjustable in length and inclination, and a special soft-rubber padding in the forehead area. **IPL protective goggles ONTOR** provide safety for users and patients during **IPL treatments**.

- suitable for **IPL applications**
- very large field of vision

Certified according to **EN 12609-1**

- Protection level: **F-4**

Scale number	Spectral transmittance $\tau(\lambda)$			Luminous transmittance $\tau_{(V,D65)}$	Near IR transmittance τ_{NIR}
	$250 \text{ nm} \leq \lambda \leq 315 \text{ nm}$ Maximum %	$315 \text{ nm} < \lambda \leq 380 \text{ nm}$ Maximum %	$380 \text{ nm} < \lambda \leq 450 \text{ nm}$ Maximum %	$380 \text{ nm} < \lambda \leq 780 \text{ nm}$ %	$780 \text{ nm} \leq \lambda \leq 1400 \text{ nm}$ Maximum %
F4	0,1	0,4	$\tau_{(V,D65)}$	$8,5 > \tau_{(V,D65)} \geq 3,2$	50

Included in delivery: microfibre bags with PROTECT-logo, polybag, user's manual, cleaning cloth with PROTECT-logo, glasses head strap